

NSF BioPACIFIC MIP

Analytical HPLC Manual



For training and further questions contact

Zach Nett: zjnett@ucsb.edu

Last Updated October 16, 2024


Contents

1	System Overview	2
2	Check the waste container	2
3	Power on the instrument	2
4	Solvent management	4
4.1	Solvent preparation	4
4.2	Solvent and additive exchange	5
5	Analytical sample preparation	6
6	Create an analytical method	8
6.1	Pump tab	8
6.2	MS tab	10
6.3	Data Acquisition tab	12
6.4	Interface tab	13
6.5	PDA tab	13
6.6	Column Oven tab	14
6.7	Autosampler tab	14
6.8	AutoPurge tab	15
6.9	Download and Save the method	15
7	Run an analytical sample	15
7.1	Start a new batch	15
7.2	Modifying a running batch	17
7.3	Manual Shutdown	17
8	Analyzing the data	18
8.1	Postrun program for analysis	18
8.1.1	Viewing operation tips	19

8.1.2	View MS chromatogram and spectra	19
8.1.3	View the PDA data	22
8.1.4	Modify the integration parameters	23
8.1.5	Export the data	24
8.2	Data Browser	24
8.2.1	Comparing multiple files	25
8.2.2	Overlay PDA and MS data	25
8.2.3	Layout Templates for comparing multiple files	26

1 System Overview

The Shimadzu Nexera analytical HPLC system within the NSF BioPACIFIC MIP facilitates the identification and quantification of complex mixtures while also supporting method development for their purification on preparative and semi-preparative scales. The system is equipped with a 6-way switch valve, allowing for automated switching between 4 columns with distinct surface chemistry and retention properties. The column selection is controlled by the programmed method. The schematic below illustrates the overall system layout.



The analytical pump is a quaternary pump that uses four different solvents located in 1 L bottles on top of the instrument. There is a 180 μL mixer located within the pump module that blends the fluids. Analytical samples are introduced by the autosampler, separated by one of four analytical columns, and subsequently analyzed by a dedicated PDA and the LCMS-2050 MS detector.

In general, running the HPLC involves six key steps:

1. Start the software and prepare the instrument.
2. Design a method specific to the analysis requirements.
3. Prepare the sample for injection.
4. Configure and execute a batch run.
5. Clean the columns and shut down the instrument after use.
6. Analyze the obtained data and export as needed.

These steps, along with other important aspects like solvent and waste management, will be further explained in the subsequent sections.


2 Check the waste container


Before starting a run, it is essential to check the capacity in the receiving waste container. If the waste container is full, it should be transferred to the designated hazardous waste pick-up area, ensuring that the waste label is correctly marked. Subsequently, a new waste container, typically an empty 4 L solvent bottle, should be installed, accompanied by a fresh hazardous waste label.


3 Power on the instrument




Each user will clean the columns that were used in their method and shutdown the instrument before leaving. When approaching the HPLC, it will be in a sleep state and look as such:




Double-click the Lab solutions icon  on the computer desktop to open the software. Login user ID is “Admin”. No password is required. Click OK.



Start the analytical  real-time analysis program by double-clicking the icon listed under Instruments. An audible beep will be heard as the system connects to the computer.

When the  icon is double-clicked, the analytical HPLC software will open to the following screen. Click  and  to ready the instrument. It will take about ten minutes for the Desolvation Temperature to reach its operating temperature of 450 °C, at which point the MS will be ready to use. To efficiently use this time, begin to set-up a method or prepare your samples. If present, the warning **LC Pump A Wait Purge** will resolve when a batch of samples is started using the batch template and autopurge function, as described in section 7.1.




4 Solvent management

A solvent manager system monitors the solvent level gravimetrically by the weight of the bottles sitting on top of the LC system. A warning will be given when the level within one bottle falls below 300 mL, and the system will shut down entirely when approximately 100 mL of that solvent remains in order to protect the columns.


4.1 Solvent preparation

Dedicated bottles of solvent employing TFA or formic acid as an additive are available for use. Before starting a run, ensure that adequate solvent is installed on top of the system and that the additives it contains are appropriate for your experiment. If more solvent is needed, this section provides instructions on how to prepare additional solvent containing 0.1% additives (formic acid or TFA) as required. LCMS grade formic acid is stored in the refrigerator, while LCMS grade TFA can be found in the same tub as the prep solvents on the floor. For transferring TFA or formic acid, utilize a disposable 1 mL syringe, as illustrated below:



To prepare additional analytical solvent, follow these steps:

1. Ensure the pump is turned off.
2. Remove the bottle to be exchanged or refilled, allowing the tubing and filter attachment to hang gently off the front side of the instrument (refer to the accompanying photos for reference).
3. Take a cap, seal the bottle selected for use, and place it on the large platform balance for taring.
4. LCMS grade solvent is stored in the blue and yellow flammable cabinet. If an open 4 L bottle is available, use it to fill the 1 L HPLC bottle with the desired solvent (e.g., water or acetonitrile) in a fume hood. If a new bottle must be opened, please label it appropriately.
5. Reseal the 1 L bottle and determine the mass of the added solvent.
6. For a 0.1% solution in water, divide the measured mass by 1000 to calculate the volume (in mL) of TFA (or formic acid) that should be added to the bottle. For a 0.1% solution in acetonitrile, divide the measured mass by 786 to calculate the milliliters of TFA (or formic acid) that should be added to the acetonitrile bottle.
7. If this bottle is the same composition as that which was removed, immediately re-install the bottle on top of the instrument. Otherwise, proceed to the next section to manually purge the pump while installing the new bottle.




4.2 Solvent and additive exchange

It is recommended to perform a manual purge of the analytical pump when switching between solvents containing formic acid, TFA, or other volatile additives. It is important to note that non-volatile additives such as phosphates or other salts should never be used on this system, as they


will damage the MS.

To purge the solvent lines, follow these steps:

1. Open the door of the analytical pump.
2. Turn the handle of the black valve counterclockwise to open the purge valve.



3. Press the purge button, type the number of the solvent line you would like to purge (e.g., 1 = Solvent A, 2 = Solvent B, 3 = Solvent C, 4 = Solvent D), and press Enter.



4. Install the new bottle of solvent on top of the system while the purge is in progress. This ensures that the new bottle remains uncontaminated by any remnants of the previous solvent in the lines.
5. Monitor the progression of air bubbles introduced during the transfer of the inlet tubing to the new bottle. Wait until all air bubbles have been eliminated after passing through the purge valve exit.
6. Repeat for each additional solvent line as needed. Once the purging is complete, remember to close the black valve finger-tight. *Do not close too tightly as it may cause damage to the pump.*


5 Analytical sample preparation

Always filter samples using the provided disposable 0.45 μm syringe filters. Both hydrophilic PVDF and hydrophobic PTFE membrane materials are suitable for filtering aqueous solutions, but sometimes more hydrophobic samples (peptoids/peptides) can adsorb to the PTFE surface when poorly solvated in water; in which case, the PVDF membranes are recommended.


Vials of the mobile phase solvents with additives are provided near the instrument to prepare samples for injection. Samples should be dissolved in 100% aqueous solution. If the sample is insoluble in 100% water, add organic modifier (acetonitrile), but keep it to a minimum to prevent premature elution of the sample from the column.

Sample concentrations of 0.001 to 0.01 mg/mL are generally adequate for injection volumes of 10–20 μL . However, ~ 0.1 mg/mL has been necessary for large peptides/peptoids (>50-mer). If you are not obtaining enough signal from your sample at a given concentration, you can opt to increase the injection volume up to 100 μL . As a rough guideline, the max sample capacity for a 4.6 mm diameter analytical column is 0.1–0.2 mg. However, it's important to note that the exact load capacity can vary depending on the nature of the sample and the specific separation conditions being used. **Avoid overloading.** Too much sample should not be injected. If peak shape changes by diluting the sample 10-fold or by decreasing the injection volume to one fifth, there is a possibility of sample overloading.

Overloading is also indicated by the fronting of peaks as illustrated:



Open the autosampler tray and place the prepared samples inside, noting the tray number and labelled position number within the tray.




6 Create an analytical method

Navigate to **File** → **New Method File**. Select **Analytical Template** to start developing a new analytical method. This template applies a 5-95% gradient of B (acetonitrile) over 45 min with a 3 min hold to flush the column and 2 min hold to re-stabilize at the initial condition.


The Method Editor window will appear, containing a series of tabs that represent different modules of the HPLC. It is recommended to modify the tabs in the order presented. Some tabs may not require modification.

6.1 Pump tab


The pump tab defines the gradient that will be applied during the sample run. When the tab is selected, it by default shows settings for pump LPGE, which refers to the quaternary analytical pump.




Analytical flow rates should be within the range of 0.5–1.2 mL/min. A gradient steepness of 0.5–2.0 % acetonitrile/min is recommended for peptides and peptoids. Time to reach flow protects the columns from rapid changes in pressure and should be set to 1 min. Under the gradient menu, a simple gradient from 5–95%B can be programmed for the method using the graphical entry on the right shown below. Always include an isocratic hold for 3 min at 95 %B at the end of the gradient to clean the column for the next run.



Gradients of greater complexity can be defined via Advanced gradient mode. An example of the above gradient in tabular form is shown below:




Use of the table format is advantageous for reducing the overall run time by applying gradients within a narrower range of interest. However, it is important to note that a subsequent isocratic hold at 95% B is always necessary at the end of the run to clean the columns. To ensure the column returns to its initial condition, a 2-minute hold at 5% B is included. An example of an advanced gradient table, which incorporates a reduced LC Stop time, a narrower eluent composition range of 5-50% B, and isocratic holds, is provided below:




6.2 MS tab

Important: Only use positive mode for solvents containing TFA or damage to the MS detector may occur. If you would like to use negative mode, ask for assistance in switching the system to an appropriate solvent and additive.

Ensure that the cycle time is set to 0.5 seconds in the MS tab. The data acquisition end time in the SCAN event table determines the duration over which the MS detector scans the m/z range and records data to generate a chromatogram of the total ion current. The acquisition time should be equal to the total run time, including any isocratic holds. Please note that the "Apply to All acquisition time" icon in the Data Acquisition tab (see section 6.3) no longer applies the LC Stop Time to the MS tab. Therefore, it is crucial to check that the end time in the event table corresponds to the stop time set in the Pump tab.




The entire m/z range for scanning is 10 to 2000. It is possible to acquire full m/z scans in both positive and negative modes simultaneously, except when TFA is used as the eluent additive. To enable negative mode, click on “Add Event” and modify the event settings to select a negative polarity and the full m/z range, as shown in the example below:




To delete an event, right-click on the entry in the event table, and select “Delete Event”.

If the analyte m/z is known, narrowing the mass range can increase instrument sensitivity. Alternatively, selective ion monitoring (SIM) can be used to enhance the sensitivity and selectivity of the analysis by focusing on specific ions of interest. In SIM mode, the mass spectrometer monitors and records the signal intensities of predetermined ions or m/z values of interest, disregarding other ions present in the sample. This targeted approach improves the signal-to-noise ratio and increases sensitivity for the desired analytes. However, it is important to note that compounds can possess multiple charges, and it is possible that the more charged compound may exhibit higher relative abundance.

To incorporate a SIM event, navigate to the SIM tab and make sure the “ m/z input” box is selected. Then, click on “Add Event” to create a new event. Choose the polarity for monitoring and use the downward arrow to modify the entries in the highlighted m/z field, which is displayed in red.




A pop-up window will appear to set the m/z values for the target analyte, as in the example shown below. Click OK to save the entries.



6.3 Data Acquisition tab


The data acquisition tab tells the detectors to acquire data for the duration of the run and to set the frequency of measurement for the analytical PDA detector. The LC Stop Time should equal the entire run time. As an example, for a 5-50% gradient applied over 25 min, LC Stop time should equal 30 min, since it includes time for cleaning and stabilization of the column.



Click **Apply to All acquisition time**. It is not recommended to adjust the sampling rates for the PDA in order to avoid data files becoming too large or the data being too sparse.


6.4 Interface tab

For the Interface tab, confirm that the drying gas is enabled as shown below.



6.5 PDA tab


The PDA is the analytical UV-Vis detector that is capable of monitoring all wavelengths (190–800 nm) simultaneously. The example below shortens the range of interest to 190–400 nm to reduce the file size and employs the deuterium (D2) lamp only. The tungsten (W) lamp or the D2 & W lamps can be used to monitor wavelengths of 371–800 nm or 190–800 nm, respectively.



6.6 Column Oven tab

Set the desired column temperature, bearing in mind that the maximum allowed temperature is 60 °C. Select the appropriate analytical column (1, 2, 3 or 4) for both valves. The column designations are as follows:

1. XBridge BEH C18 Column, 130Å, 5 µm, (4.6 mm x 150 mm)
2. XBridge BEH C18 Column, 300Å, 5 µm, (4.6 mm x 150 mm)
3. XBridge Protein BEH C4 Column, 300Å, 3.5 µm, (4.6 mm x 100 mm)
4. SeQuant ZIC-pHILIC Polymeric Column, 100Å, 5 µm, (4.6 mm x 150 mm)




For reverse-phase chromatography of small molecules and peptides/peptoids with fewer than 15 repeat units, Column 1 is recommended. Column 2 is more suitable for larger peptoids/peptides, while Column 3 is designed for highly hydrophobic materials with excessively long retention times on Column 2.

Note use of methanol in lieu of acetonitrile may increase the retention time of more hydrophilic peptoids/peptides. However, keep an eye on the increase in back-pressure when using methanol and use this solvent sparingly.

Column 4 is dedicated to HILIC (Hydrophilic Interaction Liquid Chromatography) chromatography, a technique suitable for polar and hydrophilic compounds that are challenging to retain and separate using traditional reversed-phase media. Suitable buffers for HILIC include ammonium acetate, ammonium formate, formic acid, and acetic acid, with a concentration range of 5-20 mM being sufficient for most solutes. For example, a typical gradient elution profile for the ZIC-pHILIC column is a linear gradient from 90–40% acetonitrile at a flow rate of 0.5 mL/min with a 10 min of pre-equilibration at the starting condition. The fraction of water in the eluent should not be varied outside the range of 3-60%, and the columns should never be operated at 100% water or organic. Due to the unique requirements of HILIC chromatography, it is strongly recommended to consult with BioPACIFIC MIP staff before attempting to use Column 4.


6.7 Autosampler tab

No adjustments are necessary for the settings on the Autosampler tab, which should be displayed as shown below:




6.8 AutoPurge tab

Make sure to enable the pump purge and autosampler purge options, as depicted below:



6.9 Download and Save the method


Once all the method tabs have been reviewed, click  to download the method to the system and exit the Method Editor. Save the method to your data folder by navigating to **File** → **Save Method File As...**

Create or locate your folder in the C:\LabSolutions\Data directory and name the method file such that the method is well-described. The filename “analytical_5-95ACN_45min_col1.lcm” has a good naming structure because it indicates an analytical method employing a 5–95% gradient of acetonitrile over 45 min using column 1.

7 Run an analytical sample

7.1 Start a new batch

To run a sample, select “Real-time Batch” from the main menu (refer to the image below).



Create a new batch file, by clicking **File** → **New Batch File**. Select **Analytical Batch Template**. The template includes two line entries for samples and three lines for shutdown methods corresponding to each reverse-phase column. Additional line entries can be added by inserting rows above the shutdown methods. Right-click on a line entry, select “Insert Row” to add a new row above. The shutdown methods automate flushing each column with acetonitrile so that the column is stored without TFA (or formic acid) after use. Delete shutdown method entries for unused columns. It’s important to note that using the position vial# -1 in the shutdown methods will result in no injection, effectively allowing for the running of a ”blank” without occupying a physical vial position on the sample tray.

The method template includes an AutoPurge option for the first sample run, which automatically purges the autosampler and pump. The Autopurge option should not be selected for subsequent runs.


Analysis	Tray Name	Vial#	Sample Name	Sample ID	Inj. Volume	Chromatogram Table	Analysis Type	Method File	Data File	AutoPurge
1					10		ILT		(Auto Filename)	<input checked="" type="checkbox"/>
2					10		ILT		(Auto Filename)	<input type="checkbox"/>
3	1	-1	Shutdown_Column_1		10		ILT	olutions\Data\Method Templates\shutdown_col1.lcm	(Auto Filename)	<input type="checkbox"/>
4	1	-1	Shutdown_Column_2		10		ILT	olutions\Data\Method Templates\shutdown_col2.lcm	(Auto Filename)	<input type="checkbox"/>
5	1	-1	Shutdown_Column_3		10		ILT	olutions\Data\Method Templates\shutdown_col3.lcm	(Auto Filename)	<input type="checkbox"/>

Specify the tray and vial numbers for the samples in the analytical autosampler. The sample name will be used as the filename for the data.


A typical starting injection volume is 10 μ L, but it can be adjusted from 1–100 μ L as needed. If no signal is observed from the analyte, increase the injection volume. Conversely, if the column is overloaded, decrease the injection volume. Ensure that the directory of the method matches the method file created earlier or the desired method to be used.

The example batch below shows how one can run two different samples (located in positions tray 1, vials 3 and 4). The first sample will be ran twice using two different methods - one specific for column 1 and one specific for column 2. The second sample is run also using the method developed for column 2. Because only columns 1 and 2 are used, the shutdown method entry for column 3 is deleted. To insert, add, or delete a row, right-click on the line entry and select the corresponding option.


Analysis	Tray Name	Vial#	Sample Name	Sample ID	Inj. Volume	Chromatogram Table	Analysis Type	Method File	Data File	AutoPurge
1	1	3	MWB_1-72_Col_1		10		ILT	ions\Data\Morgan\analytical_5-95B_45min_Col_1.lcm	(Auto Filename)	<input checked="" type="checkbox"/>
2	1	3	MWB_1-72_Col_2		10		ILT	ions\Data\Morgan\analytical_5-95B_45min_Col_2.lcm	(Auto Filename)	<input type="checkbox"/>
3	1	4	MWB_1-74_Col_2		10		ILT	olutions\Data\Morgan\230516_5-95D_45min_Col2.lcm	(Auto Filename)	<input type="checkbox"/>
4	1	-1	Shutdown_Column_1		10		ILT	olutions\Data\Method Templates\shutdown_col1.lcm	(Auto Filename)	<input type="checkbox"/>
5	1	-1	Shutdown_Column_2		10		ILT	olutions\Data\Method Templates\shutdown_col2.lcm	(Auto Filename)	<input type="checkbox"/>


To start the batch, highlight all rows, and click the green play button . Save the batch file in your folder using an informative filename such as the date. The system will proceed to purge the solvents and autosampler, run the samples, clean the columns, and automatically shutdown at the end of the batch.



Data can be monitored in real-time using the PDA, MS, or All viewing windows (see below). The LC tab shows changes in pressure of the analytical pump and can be useful for distinguishing whether PDA signal changes are attributable to the presence of sample or some fluctuation in pressure.



7.2 Modifying a running batch

A running batch can be paused to edit the Batch Table. To pause a batch, click the  (Edit table/Restart) icon. Set the row where the batch should pause and click OK.




Data acquisition will stop at the selected row. Add, insert, or delete entries in the batch table by right-clicking on the row and selecting the row action. **Do not delete entries that have already run.** Click the save  button on the toolbar. Remember to click the  (Edit table/Restart) icon to resume the batch.

7.3 Manual Shutdown

The instrument will shutdown automatically after a batch is completed. However, if the instrument needs to be shutdown manually, simply close the realtime analysis software window to exit the

software. A shutdown prompt will appear. Check the boxes as shown below and click OK.




8 Analyzing the data

Quantitation uses chromatographic data to determine the amount of a given component in a mixture. There are several types of quantitation methods commonly used. This section covers the very basics of data analysis, utilizing peak area percent obtained from an integrated PDA chromatogram as a rough estimate of the relative amounts of analytes present. This method assumes that the wavelength absorbance properties of all compounds observed are equivalent, which is likely not a valid assumption. The MS data is used only for mass confirmation of the observed analyte.

Other more in depth methods of quantitation require the use of internal and external standards and are not commonly employed in the routine analysis of newly synthesized materials. For more information on these types of analyses, refer to the manufacture manual entitled “Operators Guide (LCMS Edition)”.

8.1 Postrun program for analysis


Select “Postrun” from the LabSolutions Main window.





Drag-and-drop a data file onto the Data analysis window from the Data Explorer sub-window.


8.1.1 Viewing operation tips

▼ Tips About viewing operations

 **Hint** An area on a graph can be zoomed and displayed by dragging over it with the mouse. By right clicking on a graph, [Initialize Zoom] and [Undo Zoom] can be selected.

 **Hint** Drag the frame border to change the relative size of each view.


 **Hint** The display factor of the intensity axis can be incremented or decremented.



The screenshot shows a software interface with two main panels. The top panel is 'Chromatogram View' showing a chromatogram with a single prominent peak. The bottom panel is 'Spectrum View' showing a mass spectrum. Callout boxes provide tips: one points to a peak on the chromatogram, another points to the window border, and a third points to the intensity axis scale in the spectrum view.


8.1.2 View MS chromatogram and spectra

For analyzing MS data, click [MS Data Analysis] under the [Main] window.




To generate a spectrum, double-click the time point in the chromatogram. Click the arrows for “peak” or “scan” to move the spectrum extraction line to the next peak time or by the next preset scan unit.

1 Double click on the chromatogram.




Double-click the m/z in [spectrum view] to extract as a chromatogram. The chromatogram of that m/z is added to the [Chromatogram View].



Alternatively, to add extracted chromatograms of other m/z values to the chromatogram display, edit the MS Data view parameters sub-window under the MS Data tab as shown.



To average MS spectra over a specified time span:



The spectrum display can be improved by subtracting the background MS spectrum from the averaged spectrum:


The screenshot shows the PostRun Analysis software interface. A callout box labeled '1 Click' points to a toolbar icon. A second callout box labeled '2 Drag the mouse over the section of the chromatogram to be selected.' points to a peak in the chromatogram. A third callout box labeled 'Hint' points to the spectrum view below the chromatogram.

1 Click


2 Drag the mouse over the section of the chromatogram to be selected.



Hint This display indicates subtraction of the averaged spectrum between 3.883 and 4.433 minutes from the averaged spectrum between the retention times of 4.650 and 4.867 minutes.

8.1.3 View the PDA data



To view the analytical PDA in post-run, click “PDA data analysis” icon located under the main sub-menu on the left. The data should appear in a simplified form showing only the chromatogram view, peak table, and method view. If there are more windows shown, under the Layout menu, select the PDA layout, or drag the frame borders so that only these windows are shown for simplicity.



Select the  button to modify the Chromatogram, and then return to  mode to reflect any changes. While in edit mode, under the “Multi Chrom” tab, change the Ch1 wavelength of absorbance to the wavelength of interest (e.g., 214 nm for peptoids and peptides).

8.1.4 Modify the integration parameters

In edit mode, under the integration tab, select the channel of interest. If necessary, reduce the number of peaks included due to noise in the data by increasing the [Slope] or [Min. Area/Height] by a factor of 10.

Multi Chrom Integration Identification Quantitative Compound Group Performance UV

Channel: Ch1 214nm

Algorithm: Chromatopac Copy to All Channels

Width: 5 sec

Slope: 1000 uV/min

Drift: 0 uV/min


T. DBL: 1000 min

Min. Area/Height: 1000 counts


Calculated by: Area Height


Auto Max. Peak 6 counts Relative to Main Peak 1 %

Register Spectrum to Table

 Program

Noise/Drift Calculation... Advanced...

To prevent integration of extraneous peaks at early or late times in the run, click , and fill-out the table to turn off integration for specific regions:



These steps are generally sufficient for processing. Press the help button to learn more about peak integration commands if desired.

8.1.5 Export the data

To export the peak tables, chromatograms (PDA and UV), MS Chromatogram TIC, and MS Spectrum, go to **File** → **Export Data** → **Export Data (ASCII Conversion)**. Select the Output file option and check the boxes for the desired output options (peak table, chromatogram, etc.) to save the data as a tab-delimited .txt.


No MS spectrum will be exported if it has not been registered. To register a spectrum click




“register spectrum to process table” icon located in the top menu bar.

8.2 Data Browser



Open Data Browser to display chromatograms and spectra from multiple data files for comparison or to overlay data from different detectors.



8.2.1 Comparing multiple files





Within Data Browser, data files are opened by drag-and-drop from the data explorer window to the display cells on the right. To change the type of data displayed in each cell, right-click the cell and select “Change data type”.

Use the layout property icon  to change the row x column arrangement of cells. To add a row or column to an existing layout, use the  icons on the top panel.

Cells are linked by their cell number shown in the upper right-hand-corner of each cell. For example, cells labeled with “1” will correspond to the same opened file, and cells labeled “2” will correspond to a second open file. Click on the cell number to change it. Drag-and-drop files to the different sets of cell numbers to view multiple files simultaneously.

8.2.2 Overlay PDA and MS data

This section describes how to use a layout template to overlay PDA and MS data for a single analytical data file:

1. Go to the menu **Layout** → **Open Layout File**.
2. Select the layout **Analytical.lyt** from the Browser Layout Templates folder located on the Desktop. Drag and drop a file from the data explorer window to the display cells on the right.
3. Click the Data Processing Parameters icon  in the data browser menu.
4. On the “Multi Chrom” tab, enter the absorbance wavelength of interest and uncheck “Display Extracted Chromatogram”. Click OK.
5. Double-clicking within the PDA chromatogram will generate the mass spectrum for that time point. MS spectra can also be averaged or background subtracted using the same  icons described in section 8.1.2.

6. The auto-integration parameters can be adjusted (if not previously performed and saved in post-run analysis), in the Integration tab of the Data Processing parameters window.
7. To overlay the PDA and MS data within the same window. Right-click on the MS chromatogram and open “Display settings”. Navigate to the “LC Settings” tab and check the Disp. box. Click OK. The MS and PDA chromatograms are overlaid.

8.2.3 Layout Templates for comparing multiple files

Three additional layout templates have been pre-configured to assist in comparing of multiple data files:

1. **Analytical_Comparison.lyt** compares the chromatograms and spectra of two analytical data files.
2. **Prep_Comparison.lyt** compares the chromatograms and spectra two prep data files.
3. **Analytical_Prep.lyt** compares the chromatograms and spectra of one analytical and one prep data file.

Click on the menu **Layout** → **Open Layout File**. Select the layout from the Browser Layout Templates folder located on the Desktop. Drag and drop files from the data explorer window to separate cell windows labeled “1” and “2”.