


Shimadzu Nexera Preparative HPLC Manual



For training and further questions contact

Zach Nett: zjnett@ucsb.edu

Last Updated October 16, 2024


Contents

1	System Overview	2
2	Check the waste container	2
3	Solvent preparation	3
4	Start the software and power on the instrument	3
5	Purge the prep pumps	5
6	Prepare the samples	5
7	Purge the prep injector	6
8	General considerations in the use of prep HPLC	7
9	Create a prep method	7
9.1	Pump tab	8
9.2	MS tab	9
9.3	Data Acquisition tab	10
9.4	Interface, Analog Output, LC Time Prog., Controller and AutoPurge tabs	10
9.5	Detector A tab	11
9.6	Subcontroller A tab	11
9.7	Injector tab	12
9.8	Fraction Collector tab	12
9.9	Download and Save the method	16
10	Run a prep sample	16
10.1	Set up a batch	17
10.2	Modifying a running batch	19
10.3	Early termination of a run	19
10.4	Clean the columns and shut down the instrument	20

11 Analyzing the data	22
11.1 Postrun program for analysis	22
11.1.1 Viewing operation tips	23
11.1.2 View MS chromatogram and spectra	23
11.1.3 View the UV data	26
11.1.4 Modify the integration parameters	27
11.1.5 Export the data	27
11.2 Data Browser	28
11.2.1 Comparing multiple files	28
11.2.2 Overlay UV and MS data	29
11.2.3 Layout Templates for comparing multiple files	29

1 System Overview

The Shimadzu Nexera preparative HPLC system at BioPACIFIC MIP enhances the conventional UV-based fractionation process of prep chromatography by incorporating an inline single quadrupole MS detector (LCMS-2020). This integration enables direct on-line mass spectrometry, ensuring the collection of only the fractions containing the target compound and potentially eliminating the need for subsequent offline MS analysis. The system is equipped with a 6-way switch valve, allowing for automated switching between 3 semi-prep columns (19 mm ID) and 2 preparative columns (50 mm ID) that can be installed simultaneously. Column selection is conveniently controlled within the programmed method. The schematic below illustrates the overall system layout.



The preparative HPLC system operates with each prep pump supplying a single solvent to a 15 mL inline mixer, ensuring homogenization of the mobile phase. Each solvent is sourced from a 4 L solvent bottle, situated below the instrument on the floor. Samples are introduced through the injector and undergo separation by passing through one of the prep columns located in the column closet. Eluting peaks are analyzed using a dedicated UV-vis detector (referred to as Detector A). A small portion of the preparative stream is diverted to the MS detector for analysis through an inline splitter. Fraction collection can be performed manually or triggered automatically by the UV-vis detector.

In general, running the HPLC involves six key steps:

1. Start the software and prepare the instrument.
2. Design a separation method for the compound of interest.
3. Prepare the sample for injection.
4. Configure and execute a batch run.
5. Clean the columns and shut down the instrument after use.
6. Analyze the obtained data and export as needed.

These steps, along with other important aspects like solvent and waste management, will be further explained in the subsequent sections.

2 Check the waste container


Always check the solvent level in the receiving waste container prior to beginning a run. When the waste container is full, transfer the container to the hazardous waste pick-up area and check that the waste label is properly marked. Install a new waste container with a hazardous waste label attached. Depending on availability, three different types of waste containers are able to be used with the instrument.



In the left photo, tubing is connected via a hose barb adapter for the container shown. The middle photo shows a reusable white 9 L container provided by EH&S that should only be used for HPLC eluent waste. The right photo shows how to connect a recycled 4 L solvent bottle as a waste reservoir. Spare caps, with drilled holes for inserting the waste line, are located in the secondary containment tub.

3 Solvent preparation

This section describes how to prepare additional solvent containing 0.1% additives (formic acid or TFA) as needed. LCMS grade formic acid is stored in the refrigerator. LCMS grade TFA is stored in the same tub as the prep solvents on the floor to the left of the instrument. To transfer TFA or formic acid, use a disposable 1 mL syringe as depicted below.




The prep system uses 4 L bottles that are located on the floor below. Unopened 4 L bottles of LCMS grade acetonitrile or water are stored in the blue or yellow flammable cabinets located to the left of the instrument. To prepare a 0.1 vol% solution of TFA or formic acid, simply open a new bottle and add 4 mL of TFA or formic acid, using the disposable syringe. Remember to label the bottle appropriately indicating the additive used with a sharpie.


4 Start the software and power on the instrument


Each user will clean the columns that were used in their method and shutdown the instrument before leaving. When approaching the HPLC, it will be in a sleep state and look as such:






Double-click the Lab Solutions icon  on the computer desktop to open the software. Login user ID is “Admin”. No password is required. Click OK.



Start the preparative  real-time analysis program by double-clicking the icon listed under Instruments.



After double-clicking the preparative Prep LCMS icon, the real-time analysis program will open to the following screen. Click  and the four icons  to ready the instrument.



5 Purge the prep pumps

Switch the B solvent line from the 4 L pure acetonitrile bottle to the 4 L bottle of acetonitrile with 0.1% TFA (or formic acid) additive. Check that there is a sufficient amount of solvent in both the A and B bottles for your runs. If necessary, prep additional solvent (see section 3).

Purge the A and B prep pumps until any bubbles observed upon purging are removed. To purge the solvents, open the front door of the pump modules and turn the handles of the black valve counterclockwise to open the purge valves. Press the purge button.



Watch for any bubbles in the lines and let these pass through the piston heads of the pumps for removal. Press the purge button again to stop purging. Remember to close the black valve finger-tight. *Do not close too tightly or you may damage the pump.*


6 Prepare the samples

Always filter samples using the provided disposable 0.45 μm PVDF or PTFE syringe filters. Jars of the mobile phase solvents with TFA and formic acid additives are provided near

the instrument to prepare samples for injection. Samples should be dissolved in 100% aqueous solution. If the sample is insoluble in 100% water, add organic modifier (acetonitrile), but keep it to a minimum. For more hydrophobic samples that are poorly solvated in water rich solutions, the hydrophilic PVDF membranes are recommended.


Prep injection volumes can be varied between 100–2000 μL . Sample concentrations and injection volumes should be tested to determine the max sample load where adequate fraction purity is still achieved. It is best practice to maximize sample concentration rather than opting to inject a larger volume of a more dilute sample. As one example to note, peptoids have been purified with a sample load of 60 mg (e.g., 1 mL injection of a 60 mg/mL solution) on a 19 mm ID column.

Place the filtered sample in an open vial onto the prep injector tray, noting the tray number and labelled position number of the vial.




7 Purge the prep injector

Return to the software and expand the control window of the system by clicking the icon with an up arrow as indicated in the figure below.



To purge the injector, simply click on the designated icon shown below. The purging process will be completed in approximately 2 minutes.



8 General considerations in the use of prep HPLC

While analytical chromatography focuses on developing chromatography conditions that optimize resolution vs time, prep chromatography focuses on developing conditions that offer the most practical balance of yield, purity, and throughput (i.e., how much material can be purified in each injection cycle).

In analytical chromatography, small quantities of samples are loaded onto the column to ensure that the amount of sample loaded does not affect the resolution. The amount of sample that can be loaded without band broadening occurring depends on the size of the column. In contrast, in prep purification it is common practice to exceed the sample capacity and “overload” the column in order to increase yield and throughput. Prep chromatogram peaks are therefore typically broadened due to column overload, and the resulting chromatogram is not as well resolved as analytical chromatograms.

In prep chromatography, several fractions collected across the peak of interest are analyzed for impurities. Based on the analytical results, fractions with little or no impurities are pooled together. Fractions near the beginning and end of a peak which may contain significant impurities are discarded. Choosing which fractions to pool and which to discard is a matter of balancing purity and yield.

Prep chromatography uses large volumes of solvent and requires significant time. Therefore a good workflow would use analytical chromatography to first establish elution conditions that give the best separation. Fractions collected across the peak during a prep run will be re-analyzed by the analytical method to determine maximum yield and purity for pooling. Fractions are usually dilute enough to re-analyze on the analytical side without further dilution being necessary.


9 Create a prep method

Navigate to **File** → **New Method File**. Select **Prep Method Template**, and the Method Editor window will appear, containing a series of tabs that represent the modules of instrumentation. It is recommended to modify the tabs in the order presented. Some tabs may not require modification.

This template method applies a 5-95% B (acetonitrile) gradient over 45 min with a 3 min hold to flush the column and 2 min hold to re-stabilize at the initial condition.


9.1 Pump tab

The pump tab defines the gradient that will be applied during the sample run. It will by default show settings for pump B.GE1, which refers to the two prep pumps.




Flow rates ranging from 10 to 20 mL/min are suitable for the installed semi-prep columns. To maintain a similar separation quality to analytical chromatography, the prep flow rate is scaled to approximately 17 times the analytical flow rate, corresponding to the increase in column cross-sectional area. Commonly used flow rate combinations include 0.6 mL/min for analytical and 10 mL/min for prep, or 0.8 mL/min for analytical and 15 mL/min for prep. Running at 15 mL/min instead of 10 mL/min may enhance peak shape.

It is recommended to begin with an initial gradient steepness of 2.0% acetonitrile/minute. Slower gradients can be employed to enhance resolution or purity based on specific requirements. The gradient can be modified in Simple or Advanced gradient mode. Advanced gradient mode (see below) displays the programmed gradient in table form. The example below narrows the gradient to a smaller region (15–55% B) where the targeted compound elutes, thereby limiting solvent waste and reducing the run time. The required holds at 95% and 15% B for 3 and 2 min, respectively, in order to ensure the column is clean and ready for the next injection. To reduce the compositional range of the gradient, modify the table entries and reduce the “Stop time” field entry in order to shorten the total run time.




9.2 MS tab

Important: Only use positive mode for solvents containing TFA or damage to the MS detector may occur. In the MS tab, the acquisition time specifies over what time the mass spec will scan the m/z range and record data, generating a chromatogram of the total ion current (TIC). The acquisition time should equal the length of time required for the entire run, which includes time for the gradient plus 5 additional minutes to wash the column and return to the initial condition.




In the above example, the mass spectrometer is set to scan over a pre-defined mass range of 10-2000 m/z and record full mass spectra. Additional events can be added to increase detection sensitivity for a specified m/z, such as selective ion monitoring (SIM). To add a SIM event, right click on Event1Scan(+) located within the left panel, and click “Event Add”. Change the acquisition mode from scan to SIM and specify the m/z values of interest in Ch1, Ch2, etc.



9.3 Data Acquisition tab


The data acquisition tab tells the detectors over what time to acquire data and sets the frequency of measurement for the UV detector. The LC Stop Time should equal the time it takes for the run. For a 15-55% gradient applied over 20 min, LC Stop time should equal 25 min, since it includes time for cleaning and stabilization of the column.



Click **Apply to All acquisition time** so this setting will be carried over and applied to all other tabs (including the MS tab described above). It is not recommended to adjust the sampling rates.

9.4 Interface, Analog Output, LC Time Prog., Controller and AutoPurge tabs


For the interface tab, confirm that the image below matches the current settings in the method.



No changes should be necessary on the analog output, LC Time Prog., Controller, or AutoPurge tabs.

9.5 Detector A tab


Detector A is the prep UV Detector. Choose the D2 lamp for monitoring wavelengths 190–370 nm or the W lamp for wavelengths 371–1000 nm. For monitoring peptoids or peptides, use a Ch1 wavelength of 214 nm. Although two wavelengths (Ch1, Ch2) are monitored, only the Ch1 signal is used for triggering fraction collection.



9.6 Subcontroller A tab

In the Subcontroller A tab, select the same prep column for both valves in accordance with these designations:

1. XBridge BEH C18 OBD Column, 130Å, 5 µm, (19 mm x 150 mm)
2. XBridge BEH C18 OBD Column, 300Å, 5 µm, (19 mm x 150 mm)
3. XBridge Protein BEH C4 OBD Column, 300Å, 5 µm, (19 mm x 150 mm)




Column 1 is recommended for small molecules and peptides/peptoids with < 15 repeat units. Column 2 is recommended for larger peptoids/peptides. Column 3 is available for materials that are highly hydrophobic and elute at excessively long retention times on Column 2.

Note use of methanol in lieu of acetonitrile may increase the retention time of more hydrophilic peptoids/peptides. However, keep an eye on the increase in back-pressure when using methanol and use this solvent sparingly.

9.7 Injector tab

Ensure that the prep injector is enabled.



9.8 Fraction Collector tab

The fraction collector tab specifies the parameter settings for fraction collection, which can be based on three different triggers:

1. **Peak-based fractionation:** the fraction collector valve automatically opens/closes based the observation of a peak in the UV signal.

Method Parameter

Use Fraction Collector

Fraction Time: 25.00 min

Detector A Time Program Other

Use

Peak Detection Parameter

Slope

Use Slope

Front Slope: 8000 uV/sec

Back Slope: - 2000 uV/sec

Peak Shape: Unspecified

Level

Use Level

Level: 1000 uV

Slope Disable Level 10 %

Peak Collection Sensitivity: 3

Delay Time

Time Program

The example above shows default values for fraction collection based the UV signal, with the “use” boxes enabled (checked) for both the Fraction Collector and Detector A.

2. **Time-based fractionation:** the fraction collector valve is programmed to open and close at a specific retention time.

Method Parameter

Use Fraction Collector

Fraction Time: 25.00 min

Detector A Time Program Other

Perform Waveform Processing in the Lock Section

	Time	Command	Parameter
1	8.00	ValveOpen	
▶ 2	8.50	ValveClose	▼
3			
4			
5			

3. **Manual fractionation:** during a run, the user manually opens and closes the fraction collector valve based on the visual observation of a peak in the LC data output window. To manually open/close valves, use these equivalent sets of buttons located at the top of software menu or within the control panel, respectively:




Buttons located in control panel window below data output display. Buttons will be active during a run.




Note that there will be a time delay of ~ 11 seconds between the click and opening of the valve due to the position of the fraction collector relative to the UV detector.

The collection simulation window allows the user to confirm the appropriateness of the collection starting and stopping criteria. This is best used by opening a data file generated during a prep “test run” that involved only a small injection volume of 100 μ L. Collection simulation is an optional software feature used to tune the parameter settings for peak-based fractionation:

1  **2** Click [Add] to add target peaks for preparative work to the detector chromatogram. Peaks can be inserted by clicking to select the starting and stopping positions.

3 Click [Reflect]. Preparative parameter settings for [Front Slope], [Back Slope], and [Level] are generated automatically.


NOTE In some cases, peaks that were not specifically added are identified by selecting peak starting and stopping positions. If unwanted peaks are recognized, use Lock/Unlock peak integration commands in a time program to ensure peak integration does not occur for the unwanted peaks.




To use collection simulation to tune the parameter settings for time-based fractionation:

2 Select the time program to add and select the chromatogram range in preparative simulation results.


- Hint**
- Peak integration Lock/Unlock commands cannot be added to intervals with the same command already specified.
 - Added programs can be deleted using mouse operations.



3 Click [Reflect].
Time programs are generated automatically based on preparative parameters.
Time programs added based on preparative simulation results are applied to preparative parameter time programs.




9.9 Download and Save the method

Once all the method tabs have been reviewed, click  to download the method to the system and exit the Method Editor. Save the method to your data folder by navigating to **File** → **Save Method File As...**

Create or locate your folder in the C:\LabSolutions\Data directory and name the method file such that the method is well-described. The filename “prep_5-95ACN_45min_col1.lcm” has a good naming structure because it indicates a prep method employing a 5–95% gradient of acetonitrile over 45 min using column 1.


10 Run a prep sample

Ensure that all collection test tubes are clean, empty, and in place on the fraction collector. Check the rack status in the software, and clear the racks in the control panel of the fraction collector if necessary (see below).



10.1 Set up a batch

To run a sample, select “Real-time Batch” located under the main menu.



Create a new batch file, by clicking **File** → **New Batch File**. Select **Prep Batch Template**. Fill out the table with line entries for each injection/sample. It is recommended to include at least one blank injection (using vial position “-1”) to condition the columns prior to a sample run. To insert (or add) a row, highlight a line entry and right click. Select “Insert Row” (or “Add Row”). Input the following information for each sample/injection as defined below.


Analysis	Tray Name	Vial#	Sample Name	Sample ID	Method File	Data File	Inj. Volume	UV Wavelength	Collection Level	Fraction Vessel	Divide Volume
1	1	-1	No Injection_system warmup		od Files\Prep_30-60ACH_20min_col4.lcm	(Auto Filename)	10				
2	3	2	MWB_1-35	500UL	od Files\Prep_30-60ACH_20min_col4.lcm	(Auto Filename)	500				

These inputs will default to method values if left blank


- **Tray Name and Vial#** = position of the sample vial in prep injector tray. “-1” vial number results in no injection, i.e. “a blank sample”.



- **Sample Name** = name of generated data file.
- **Sample ID** = entry not required but is appended to the filename when used.
- **Method File** = check that this path matches the name of the method created (or other method file you would like to use).
- **UV Wavelength** = the wavelength of absorbance used to trigger collection based on peak detection in signal (otherwise uses Ch1 specified in Detector A tab of method file).
- **Level** = specifies the level in used μV for peak detection in UV detector output.
- **Fraction Vessel** = sets the starting (or ending) rack or tube to be used in fraction collection.
- **Divide Volume** = volume used for each fraction collected. Setting range is 0-25 mL.

Preparative columns are more sensitive to changes in pressure than analytical columns. Therefore, **a batch should never use two different columns in a continuous sequence of runs.** If pump flow is not stopped before switching columns, the column can be damaged. Please contact BioPACIFIC MIP staff if you need to switch between columns during a batch run, and they will help you set-up a transition method.


To start the batch, highlight all rows to be run, and click the green play button . Confirm the batch execution range and click Start.

Analysis	Tray Name	Val#	Sample Name	Sample ID	Method File	Data File	Inj. Volume	UV Wavelength	Collection Level	Fraction Vessel	Divide Volume
1		1	No Injection system warmup		od Files\Prep 30-60ACN 20min cold.lcm	(Auto Filename)	10				
2		2	MWB 1-35	500L	od Files\Prep 30-60ACN 20min cold.lcm	(Auto Filename)	500				





Save the batch file in your folder using an informative filename such as the date. The batch template is preconfigured to automatically shutdown the pump at the end of a batch in order to save solvent. If the pump is shut-off at the conclusion of a batch, click  to re-activate the pump. **Remember to stop the pumps at the end of a batch by clicking  should they fail to turn off automatically!**



During a run, the data can be monitored in the LC, MS, or All viewing windows/tabs. To change the display, right click on the window to modify the “Graph Properties” or “Display settings”.





10.2 Modifying a running batch


A running batch can be paused to edit the Batch Table. To pause a batch, click the  (Edit table/Restart) icon. Set the row where the batch should pause and click OK.



Data acquisition will stop at the selected row. Add, insert, or delete entries in the batch table by right-clicking on the row and selecting the row action. **Do not delete entries that have already ran.** Click the save  button on the toolbar. Click the  (Edit table/Restart) icon to resume the batch.

10.3 Early termination of a run

Data acquisition can be stopped at any time by selecting either stop icon  or . Check both boxes and click OK to end the run immediately. **Remember to flush the column with 95–100% B to remove any residual material remaining on the column before continuing with subsequent injections.** This can be done in the control panel by changing the Pump System A B. Conc:



Click apply and allow the pump to continue to run for ~3 min. Remember to turn off the pumps after flushing the columns.

10.4 Clean the columns and shut down the instrument

Each prep column used should be flushed with 10 column volumes of pure acetonitrile and stored without TFA. To automate the cleaning procedure, so-called shutdown methods have been created to facilitate cleaning of each column and shutdown of the instrument.

To start the shutdown procedure, switch the solvent line of prep pump B back to the bottle of pure acetonitrile. Ensure there is sufficient solvent for the cleaning procedure (~400 mL per column). Add one entry line to the batch table for each column that needs to be cleaned. The example below shows an entry for cleaning prep column 4:

Analysis	Tray Name	Vial#	Sample Name	Sample ID	Method File	Data File	Inj. Volume
1	1	-1	No Injection_system warmup		Iutions\Data\Morgan\Method Files\Prep_30-60ACN_20min_col4.lcm	(Auto Filename)	10
2	3	2	MWB_1-35	500uL	Iutions\Data\Morgan\Method Files\Prep_30-60ACN_20min_col4.lcm	(Auto Filename)	500
3	3	-1	clean_column_4		C:\LabSolutions\Data\Method Templates\Prep Shutdown Col_4.lcm	(Auto Filename)	500

The tray name does not matter. The vial# must be -1 so that no injection occurs. The method file selected corresponds to the column that will be cleaned. Shutdown method files are located in the method templates folder.

Next, right-click within the batch window and navigate to settings.

Folder: C:\LabSolutions\Data

Analysis	Tray Name	Vial#	Sample Name	Sample ID	Method File
1	1	-1	No Injection_system wamup		lutions\Data\Morgan\Method Files\Prep_30-60ACN_20min_col4.lcm
2	3	2	MWB_1-35	500uL	lutions\Data\Morgan\Method Files\Prep_30-60ACN_20min_col4.lcm
3	3	-1	clean_column_4		C:\LabSolutions\Data\Method Templates\Prep Shutdown Col_4.lcm

Within the shutdown tab, fill-out the boxes as shown below and click OK.

Shutdown

Shutdown

Shutdown Method File:
 C:\LabSolutions\Data\Method Templates\Prep Shutdown Col_5.lcm

Cool down Time: 0 min
 (After the shutdown method is downloaded, the cooldown is performed during the cool down time. After the cool down time passes, the shutdown is performed.)

LC

Degassing Unit Off
 LC Pump Off
 LC Detector Off
 Power Off after shutdown

MS

IG Off
 Nebulizing Gas Off
 DL Heater Off
 Heat Block Off
 Dry Gas Off

OK Cancel Help

Save and start the batch for the shutdown method. The system will now proceed to clean the column (~ 20 min) and then shutdown automatically.


11 Analyzing the data

Quantitation uses chromatographic data to determine the amount of a given component in a mixture. There are several types of quantitation methods commonly used. This section covers the very basics of data analysis, utilizing peak area percent obtained from an integrated UV or PDA chromatogram as a rough estimate of the relative amounts of analytes present. This method assumes that the wavelength absorbance properties of all compounds observed are equivalent, which is likely not a valid assumption. The MS data is used only for mass confirmation of the observed analyte.

Other more in depth methods of quantitation require the use of internal and external standards and are not commonly employed in the routine analysis of newly synthesized materials. For more information on these types of analyses, refer to the manufacture manual entitled “Operators Guide (LCMS Edition)”.

11.1 Postrun program for analysis


Select “Postrun” from the LabSolutions Main window.





Drag-and-drop a data file onto the Data analysis window from the Data Explorer sub-window.


11.1.1 Viewing operation tips

▼ Tips About viewing operations

 **Hint** An area on a graph can be zoomed and displayed by dragging over it with the mouse. By right clicking on a graph, [Initialize Zoom] and [Undo Zoom] can be selected.

 **Hint** Drag the frame border to change the relative size of each view.


 **Hint** The display factor of the intensity axis can be incremented or decremented.



The screenshot shows two windows: 'Chromatogram View' and 'Spectrum View'. The chromatogram shows a single sharp peak at approximately 4.767 minutes. The spectrum view shows a plot of intensity versus m/z. A 'Spectrum Type' dialog box is open, showing options for 'Centroid', 'Profile', 'Absolute Inten.', and 'Relative Inten.'. A 'Stack' button is visible in the bottom right corner of the spectrum view.


11.1.2 View MS chromatogram and spectra

For analyzing MS data, click [MS Data Analysis] under the [Main] window.




To generate a spectrum, double-click the time point in the chromatogram. Click the arrows for “peak” or “scan” to move the spectrum extraction line to the next peak time or by the next preset scan unit.

1 Double click on the chromatogram.




Double-click the m/z in [spectrum view] to extract as a chromatogram. The chromatogram of that m/z is added to the [Chromatogram View].




Alternatively, to add extracted chromatograms of other m/z values to the chromatogram display, edit the MS Data view parameters sub-window under the MS Data tab as shown.




To average MS spectra over a specified time span:



The spectrum display can be improved by subtracting the background MS spectrum from the averaged spectrum:




1 Click 




2 Drag the mouse over the section of the chromatogram to be selected.

Hint This display indicates subtraction of the averaged spectrum between 3.883 and 4.433 minutes from the averaged spectrum between the retention times of 4.650 and 4.867 minutes.

11.1.3 View the UV data

To view the prep UV chromatogram in post-run, click “Data Analysis” icon  located in the main sub-menu on the left. Select the  Edit button to modify the Chromatogram, and then return to  View mode to reflect any changes.



Peak#	Ret. Time	Area	Height	Area%
1	11.213	139468	7500	0.480
2	12.515	126645	8021	0.436
3	13.255	273171	15461	0.940
4	13.624	128319	9571	0.437
5	13.876	195681	14121	0.673
6	14.421	434239	20801	1.493
7	14.816	504096	22546	1.736
8	15.228	509090	33432	1.751
9	15.705	5935859	286413	20.415
10	16.246	541428	30765	1.862
11	16.776	489422	24340	1.652
12	17.293	338411	18161	1.164
13	17.683	306482	17532	1.054
14	18.039	30562	15266	1.049
15	18.346	296439	19482	1.020
16	18.599	315247	20970	1.084
17	18.792	251651	19543	0.866
18	19.288	1249616	90011	4.298

11.1.4 Modify the integration parameters

In edit mode, under the integration tab, select the channel of interest. If necessary, reduce the number of peaks included due to noise in the data by increasing the [Slope] or [Min. Area/Height] by a factor of 10.

Multi Chrom Integration Identification Quantitative Compound Group Performance UV

Channel: Ch1 214nm

Algorithm: Chromatopac Copy to All Channels

Width: 5 sec

Slope: 1000 uV/min

Drift: 0 uV/min


T. DBL: 1000 min


Min. Area/Height: 1000 counts

Calculated by: Area Height

Auto Max. Peak 6 counts Relative to Main Peak 1 %

Register Spectrum to Table

To prevent integration of extraneous peaks at early or late times in the run, click , and fill-out the table to turn off integration for specific regions:




These steps are generally sufficient for processing. Press the help button to learn more about peak integration commands if desired.

11.1.5 Export the data


To export the peak tables, UV chromatogram, MS Chromatogram TIC, and MS Spectrum, go to **File** → **Export Data** → **Export Data (ASCII Conversion)**. Select the Output file option and

check the boxes for the desired output options (peak table, chromatogram, etc.) to save the data as a tab-delimited .txt.


Note that no MS spectrum will be exported if it has not been registered. To register a spectrum click “register spectrum to process table” icon  located in the top menu bar.

11.2 Data Browser



Open Data Browser to display chromatograms and spectra from multiple data files for comparison or to overlay data from different detectors.



11.2.1 Comparing multiple files




Within Data Browser, data files are opened by drag-and-drop from the data explorer window to the display cells on the right. To change the type of data displayed in each cell, right-click the cell and select “Change data type”.


Use the layout property icon  to change the row x column arrangement of cells. To add a row or column to an existing layout, use the  icons on the top panel.

Cells are linked by their cell number shown in the upper right-hand-corner of each cell. For example, cells labeled with “1” will correspond to the same opened file, and cells labeled “2” will correspond to a second open file. Click on the cell number to change it. Drag-and-drop files to the different sets of cell numbers to view multiple files simultaneously.

11.2.2 Overlay UV and MS data

This section describes how to use a layout template to overlay UV and MS data for a single prep data file:

1. Go to the menu **Layout** → **Open Layout File**.
2. Select the layout **Prep.lyt** from the Browser Layout Templates folder located on the Desktop. Drag and drop a file from the data explorer window to the display cells on the right.
3. Double-clicking within the PDA chromatogram will generate the mass spectrum for that time point. MS spectra can also be averaged or background subtracted using the same  icons described in section 11.1.2.
4. The auto-integration parameters can be adjusted (if not done previously performed and saved in post-run analysis), using the Integration tab of the Data Processing Parameters window.



The  icon is found in the data browser menu on the left-hand side.

5. To overlay the UV and MS data within the same window. Right-click on the MS chromatogram and open “Display settings”. Navigate to the “LC Settings” tab and check the Disp. box. Click OK. The MS and UV chromatograms are overlaid.

11.2.3 Layout Templates for comparing multiple files

Three additional layout templates have been pre-configured to assist in comparing of multiple data files:

1. **Analytical_Comparison.lyt** compares the chromatograms and spectra of two analytical data files.
2. **Prep_Comparison.lyt** compares the chromatograms and spectra two prep data files.
3. **Analytical_Prep.lyt** compares the chromatograms and spectra of one analytical and one prep data file.

Click on the menu **Layout** → **Open Layout File**. Select the layout from the Browser Layout Templates folder located on the Desktop. Drag and drop files from the data explorer window to separate cell windows labeled “1” and “2”.