

# Diamond Saw

Diamond Saw



Tool Type: "Cutting Tool"

Location: "Microfluidics Lab"

| Supervisor            | Tool Lead |
|-----------------------|-----------|
| David Bothman         |           |
| (805) 893-4125        |           |
| bothman@cnsi.ucsb.edu |           |

Description: "Trim Saw 2"

Manufacturer: "Allied"

## About

The Trimsaw 2 is a saw designed for rapid sectioning of printed circuit boards, ceramic substrates and other abrasive non metals. Contains a built in reservoir which can be used to feed lubrication directly to the blade.

## Training

## Detailed Specifications

Motor power: .25 HP Speed range: 500-3000 RPM by 100 RPM increments Cutting platform: 6.75" x 12.5" Blade size: 4"-6"

## Safety Concerns

- Keep hands clear of blade when machine is on
  - Use only approved blades and do not cut metals with sectioning blades
  - Use coolant or lubrication when making cuts
  - Use fixtures and guides whenever possible to reduce risk of slipping into the blade
- 

## Reference Documentation

---

From:

<https://bpm-wiki.cnsi.ucsb.edu/> - NSF BioPACIFIC MIP Wiki

Permanent link:

[https://bpm-wiki.cnsi.ucsb.edu/doku.php?id=diamond\\_saw&rev=1602607814](https://bpm-wiki.cnsi.ucsb.edu/doku.php?id=diamond_saw&rev=1602607814)

Last update: **2020/10/13 16:50**

