2025/11/06 03:52 1/1 EinScan-SP Training SOP

# **EinScan-SP Training SOP**

### **Tool Name and Location**

The SHINING3D EinScan-SP 3D Scanner is located in Fume Hood #1 in Elings 2448.

# **Safety**

#### **Machine Limitations:**

• The turntable has a maximum weight limit of 5kg (~11lbs). Exceeding this weight could damage the machine and its internal components.

#### Environmental:

• Keep the scanner in a well ventilated area. Keep flammables, corrosive agents, and liquids away from the scanner. Failing to do so will damage the machine and its internal components.

## **Physical:**

• The projector emits a bright light. Avoid looking directly into the light source to prevent eye injury.

#### From:

https://bpm-wiki.cnsi.ucsb.edu/ - NSF BioPACIFIC MIP Wiki

Permanent link:

https://bpm-wiki.cnsi.ucsb.edu/doku.php?id=einscantrainginsop&rev=1597690131

Last update: **2020/08/17 18:48** 

