

PlayGround

UCLA Living BioFoundry

Fluent/Momentum Protocol

*Important:

1. Check inventory in momentum - make sure physical matches w/software
2. Check worktable in fluent - make sure physical deck mimics what software says
3. Check fluent control variables - make sure variable values are all "1" (i.e if starting at beginning of protocol)

From:
<https://bpm-wiki.cnsi.ucsb.edu/> - NSF BioPACIFIC MIP Wiki

Permanent link:
<https://bpm-wiki.cnsi.ucsb.edu/doku.php?id=playground:playground&rev=1693614684>

Last update: **2023/09/02 00:31**

