UCLA Living BioFoundry

PlayGround

Fluent/Momentum Protocol

*Important:

1. Check inventory in momentum - make sure physical matches w/software

2. Check worktable in fluent - make sure physical deck mimics what software displays

3. Check fluent control variables - make sure variable values are all "1" (i.e if starting at beginning of protocol)

From: https://bpm-wiki.cnsi.ucsb.edu/ - NSF BioPACIFIC MIP Wiki

Permanent link: https://bpm-wiki.cnsi.ucsb.edu/doku.php?id=playground:playground&rev=1693614740

Last update: 2023/09/02 00:32

