

ELN: Well Plate Operations -- Annotating Individual Wells

Biological well plates in the ELN can be annotated on a well-by-well basis by clicking on a well in the well plate diagram, which subsequently opens an individual well annotation window to the right of the plate diagram, as shown in Fig 1.



Fig 1. Well plate diagram with individual well annotation window shown on the right side.

From:

<https://bpm-wiki.cnsi.ucsb.edu/dokuwiki/> - NSF BioPACIFIC MIP Wiki

Permanent link:

<https://bpm-wiki.cnsi.ucsb.edu/dokuwiki/doku.php?id=eln-guide-basic-individual-well-annotations&rev=1744049936>

Last update: **2025/04/07 18:18**

