

## ELN: Well Plate Operations -- Annotating Individual Wells

Biological well plates in the ELN can be annotated on a well-by-well basis by clicking on a well in the well plate diagram, which subsequently opens an individual well annotation window to the right of the plate diagram, as shown in Fig 1.



*Fig 1. Well plate diagram with individual well annotation window shown on the right side.*

From:

<https://bpm-wiki.cnsi.ucsb.edu/dokuwiki/> - NSF BioPACIFIC MIP Wiki

Permanent link:

<https://bpm-wiki.cnsi.ucsb.edu/dokuwiki/doku.php?id=eln-guide-basic-individual-well-annotations&rev=1744049936>

Last update: **2025/04/07 18:18**

