

## About



Thermo Fisher TSQ Altis inline triple quadrupole mass spectrometer coupled with an ultra-high-performance liquid chromatograph (UHPLC/MS/MS) for separation and analysis of synthesized bio-derived monomers. The system is ideal for targeted metabolite analysis. The device can accurately measure between 5-2000  $m/z$  mass range with Active Ion Management (AIM) technology. The device also contains 6 channel high-pressure solvent blending, allowing a variety of solvents to run through various columns depending on the compound.

The system offers high performance integration, high-throughput capabilities, enhanced sensitivity and speed, and application versatility. With the system, researchers can acquire high-throughput and high-sensitivity targeted quantification while still maintaining robust performance across a wide variety of applications

## Reference Documentation

tsq\_hardware\_manual.pdf

trainingmaterial.pptx

### Removing and Cleaning Ion Transfer Tube:

removing\_and\_cleaning\_ion\_transfer\_tube.pdf

From:

<https://bpm-wiki.cnsi.ucsb.edu/> - NSF BioPACIFIC MIP Wiki

Permanent link:

<https://bpm-wiki.cnsi.ucsb.edu/doku.php?id=lbf-tsqaltis>

Last update: 2025/06/27 18:11



