

PlayGround

UCLA Living BioFoundry

Fluent/Momentum Protocol

*Important:

1. Check inventory in momentum - make sure physical matches w/software
2. Check worktable in fluent - make sure physical deck mimics what software displays
3. Check fluent control variables - make sure variable values are all "1" (i.e if starting at beginning of protocol)

From:

<https://bpm-wiki.cnsi.ucsb.edu/dokuwiki/> - NSF BioPACIFIC MIP Wiki

Permanent link:

<https://bpm-wiki.cnsi.ucsb.edu/dokuwiki/doku.php?id=playground:playground&rev=1693614740>

Last update: **2023/09/02 00:32**

