Vapourtec Flow Chemistry System

RS-400 with Fraction Collector png Tool Type: Scalable synthesis

Manufacturer: Vapourtec
Location: Elings Hall 2411

Principal ScientistMorgan Bates

morganbates@ucsb.edu

About

The Vapourtec flow chemistry R-series system is a flexible, precise, and automatable platform designed for efficient synthesis, reaction optimization, and scale-up. The user-friendly walk-up interface allows all users to leverage the advantages of flow chemistry, including precise temperature control, shortened reaction times, low risk in handling dangerous compounds, and scalability. The Vapourtec system at BioPACIFIC features a range of reactors offering diverse conditions for thermal, photochemical, electrochemical, and heterogeneous reactions.

Other key features include:

- · Automated reagent addition
- Fraction collection triggered by inline analytics including UV-Vis, NMR, and FT-IR
- FlowIR FT-IR analysis from Mettler Toledo
- 60 MHz ¹H NMR flow-cell based analysis from Magritek no deuterated solvents necessary!
- UV-vis analysis with variable wavelength (190-750 nm)
- Photochemistry supported at 365, 400, and 450 nm wavelengths
- Electrochemistry with broad electrode compatibility and constant current or voltage modes
- Temperatures spanning -70 to 250 °C and wide-ranging flow rates (0.02-10 mL/min)

Reference Documentation

Electrochemical Manual R2 Operator Manual

UV-150 User Manual

update: 2024/10/15 vapourtec_flow_chemistry_system https://bpm-wiki.cnsi.ucsb.edu/doku.php?id=vapourtec_flow_chemistry_system&rev=1728954002 01:00

From:

https://bpm-wiki.cnsi.ucsb.edu/ - NSF BioPACIFIC MIP Wiki

https://bpm-wiki.cnsi.ucsb.edu/doku.php?id=vapourtec_flow_chemistry_system&rev=1728954002

Last update: 2024/10/15 01:00

